

The Baxter logo is centered at the top of the slide. It consists of the word "Baxter" in a bold, italicized, blue sans-serif font. The background of the slide is a complex geometric pattern of overlapping triangles and lines in various shades of blue, gold, and light gray, creating a dynamic, abstract design.

Baxter

TMP and Filter Pressure Drop

Monitoring Circuit Pressure Trends

Presentation Overview

- Importance of therapy continuity for delivery of adequate CRRT dose
- Impact of circuit down-time on CRRT therapy continuity
- Clogging and clotting as primary causes of circuit changes
- Factors affecting CRRT circuit life
- Overview of CRRT circuit hemodynamics
- Transmembrane pressure (TMP) and filter pressure drop (ΔP filter)
- TMP and pressure drop trends as indicators of filter/circuit function
- TMP and filter pressure drop troubleshooting
- Key takeaways
- References

Therapy Continuity and Circuit Down-Time

Therapy Continuity Is Important

- One of the greatest challenges of CRRT is ensuring treatment continuity¹
- Frequent treatment interruptions can lead to significant differences between prescribed and delivered doses of CRRT, resulting in decreased therapy dosing that may have an affect on morbidity and mortality¹⁻²

1. Dunn WJ & Sriram S. Crit Care Resusc. 2014 Sep;16(3):225-31.
2. Ronco C et al. Lancet. 2000 Jul 1;356(9223):26-30.

Randomized Study of CRRT Dose and Survival²




Figure 2: Kaplan Meier estimation of survival rates in the three groups

Group 1: 20 mL h⁻¹ kg⁻¹
Group 2: 35 mL h⁻¹ kg⁻¹
Group 3: 45 mL h⁻¹ kg⁻¹

Down-Time Reduces Prescription Delivery

- Prospective study of down-time and its impact on uremic control
- 78% of down-time was due to filter clotting
- Median down-time was 3 hours, resulting in patients receiving approximately 78% of the prescribed dose
- Down-time was significantly correlated with increased plasma creatinine and urea concentrations ($p < 0.0001$)
- Down-time should be < 8 hours to maintain adequate treatment




Fig. 1 Scattergram of down-time and Δ creatinine. Univariate analyses showed significant relationship between down-time and Δ creatinine ($p < 0.0001$). $\% \Delta$ creatinine = $0.13 - 0.015 \times$ down-time; $R^2 = 0.153$




Fig. 2 Scattergram of down-time and Δ urea. Univariate analyses showed significant relationship between down-time and Δ urea ($p < 0.0001$). $\% \Delta$ urea = $0.148 - 0.02 \times$ down-time; $R^2 = 0.219$

Clogging, Clotting & Circuit Changes

- Most circuit changes are related to membrane clogging and clotting¹⁻⁶
- Frequent filter changes contribute to:
 - Incomplete dose/ prescription delivery¹⁻⁶
 - Decreased solute, fluid balance and acid-base control²⁻³
 - Increased blood loss⁶
 - Increased nursing workload⁶
 - Increased costs⁶

1. Ejaz AA et al. Nurs Crit Care. 2007 Mar-Apr;12(2):81-5.

2. Dunn WJ & Sriram S. Crit Care Resusc. 2014 Sep;16(3):225-31.

3. Fealy N et al. Crit Care Resusc. 2002 Dec;4(4):266-70.

4. Kim IB et al. Blood Purif. 2010;30(2):79-83.

5. Michel T et al. Curr Opin Crit Care. 2018 [Epub ahead of print]

6. Zhang L et al. Blood Purif. 2016;41(4):254-63.

Factors Affecting Circuit Life¹⁻²

Patient Factors	<ul style="list-style-type: none">• Pathology (sepsis, inflammation, coagulopathy)• Hematocrit• Platelet function/count
Filter Characteristics	<ul style="list-style-type: none">• Adsorptive vs. non-adsorptive• Surface area• Pore size and number
Catheter Properties & Placement	<ul style="list-style-type: none">• Lumen diameter, length, material, coating• Location
Anticoagulation	<ul style="list-style-type: none">• Anticoagulation vs. no anticoagulation• Heparin vs. regional anticoagulation
Filtration Fraction	<ul style="list-style-type: none">• higher filtration fractions associated with more clogging
Vascular Access Quality	<ul style="list-style-type: none">• Stasis, blood flow rate
Replacement Fluids	<ul style="list-style-type: none">• Pre-dilution vs. post-dilution

1. Joannidis M & Oudemans-van Straaten HM. Crit Care. 2007;11(4):218.

2. Brain M et al. BMC Nephrol. 2017 Feb 20;18(1):69.

CRRT Circuit Hemodynamics

CRRT Circuit Hemodynamics


- Understanding circuit hemodynamics will enable you to¹
 - Respond to alarms generated by CRRT devices during therapy
 - Recognize common issues like catheter dysfunction and filter clogging or clotting
- Slow or ineffective troubleshooting may result in the need for frequent circuit changes¹⁻²
- Frequent circuit changes decrease therapy efficacy and increase nursing workload and CRRT costs¹⁻²

1. Michel T et al. Curr Opin Crit Care. 2018 [Epub ahead of print]


2. Dunn WJ & Sriram S. Crit Care Resusc. 2014 Sep;16(3):225-31.

Overview: CRRT Circuit Hemodynamics

Measured Pressures




Software-Calculated Pressures



Access Pressure

Access Pressure



Pressure created by pulling blood from patient

Can be negative or positive, depending on the blood source to which the access line is connected¹

Main determinants:²

- Blood flow rate
- Catheter diameter, length, design
- Catheter tip position
- Catheter patency
- Patient hemodynamics
- Patient position and movements

1. Prismaflex Operator's Manual.

2. Michel T et al. Curr Opin Crit Care. 2018 [Epub ahead of print]

Filter Pressure


Filter Pressure

Pressure to push blood through filter

Always positive and higher than return pressure¹

Main determinants:²

- Blood flow rate
- Number of patent filter capillaries
- Filter size
- Hematocrit
- Blood viscosity
- Outflow pressure (transmitted)



1. Prismaflex Operator's Manual.

2. Michel T et al. Curr Opin Crit Care. 2018 [Epub ahead of print]


Effluent Pressure

Effluent Pressure

Can be either positive or negative, depending on whether fluid is being pushed or pulled from the blood compartment to the fluid compartment¹

Main determinants:²

- Set ultrafiltration rate
- Membrane pore permeability



1. Prismaflex Operator's Manual.

2. Michel T et al. Curr Opin Crit Care. 2018 [Epub ahead of print]

Return Pressure

Return Pressure


Pressure created by returning blood to patient

Always positive¹

Reflects resistance/pressure needed to overcome in order to return blood to patient²

Main determinants:²

- Blood flow rate
- Catheter section, length, design
- Catheter tip position
- Catheter patency
- Patient hemodynamics
- Clotting in deaeration chamber
- Patient position and movements



1. Prismaflex Operator's Manual.

2. Michel T et al. Curr Opin Crit Care. 2018 [Epub ahead of print]

Software-Calculated Pressures

Transmembrane Pressure (TMP)

Pressure gradient across the filter membrane – difference in pressure between between blood and dialysis compartments¹

Main determinants:¹


- Membrane pore permeability
- Hydrostatic/oncotic pressure gradients across filter
- Rate of replacement solutions
- Blood flow rate

Filter Pressure Drop (ΔP Filter)

Difference between Return and Filter pressures

Main determinants:¹

- Number of patent capillaries in the filter
- Patient blood viscosity



1. Michel T et al. Curr Opin Crit Care. 2018 [Epub ahead of print]

Transmembrane Pressure and Pressure Drop

Transmembrane Pressure (TMP)


$$\text{TMP} = \frac{\text{Filter Pressure} + \text{Return Pressure} - \text{Effluent Pressure}}{2}$$

Transmembrane Pressure (TMP) is the pressure exerted on the filter membrane during CRRT¹

- Reflects the pressure difference between the blood and fluid compartments of the filter

During treatment, permeability of the membrane decreases due to protein coating (clogging), causing TMP to increase

- TMP above +300 mmHg will produce the “Advisory: TMP Too High” alarm
- If the TMP increases beyond the membrane capacity that is product dependent, the “Caution: TMP Excessive alarm” occurs



1. Prismaflex Operator's Manual.

Filter Pressure Drop (ΔP Filter)


Filter Pressure Drop (ΔP Filter) = Filter Pressure — Return Pressure

ΔP Filter reflects pressure conditions in the hollow fibers

ΔP Filter is directly proportional to resistance at a given blood flow rate

Microclotting occurs in the hollow fibers during treatment, eventually leading to gross clotting and the need to change to a new set

- Clotting creates resistance to blood flow through the filter fibers and causes ΔP Filter to rise
- The amount of increase above the initial ΔP Filter contributes to the “Advisory: Filter is Clotting” alarm



Monitoring TMP and Pressure Drop Trends

TMP & ΔP Filter: Indicators of Circuit Function

TMP


- A gradual rise in TMP indicates that the filter is starting to **CLOG**
 - TMP increases as filter pores clog
- ΔP Filter
- A gradual rise in ΔP Filter indicates that the filter is starting to **CLOT**
- ΔP Filter increases as clotting creates resistance through the hollow filter fibers

ΔP Filter

- A gradual rise in ΔP Filter indicates that the filter is starting to **CLOT**
 - ΔP Filter increases as clotting creates resistance through the hollow filter fibers

What Do TMP and ΔP Filter Trends Tell You?

- Straight lines = therapy is running and filter is efficient
- Gradual increase in TMP (green line) = filter is most likely CLOGGING
- Gradual increase in ΔP Filter (white line) = filter is most likely CLOTTING
- The faster the increase, the faster the filter is degrading
- Abrupt peaks and troughs indicate therapy or blood pump has stopped
 - Access or return issues, changing a bag, etc.



Remember...

The operating points for TMP and Pressure Drop are recalculated during treatment

Software sets the initial TMP value when initial pressure operating points are established

- TMP value is subsequently reset each time the blood flow, patient fluid removal or replacement solution rates are changed, as well after self-test

Software sets the initial ΔP Filter value when the initial pressure operating points are established

- ΔP Filter value is subsequently reset each time the blood flow rate is changed

Troubleshooting: TMP and ΔP Filter

TMP and ΔP Filter Troubleshooting

	Potential Cause†	Potential Intervention†
High TMP	Membrane clogging	<ul style="list-style-type: none"> • Plan circuit replacement • Increase blood flow rate • Decrease replacement rate
	Filter clotting	<ul style="list-style-type: none"> • Plan circuit replacement • Increase blood flow rate • Increase pre-filter replacement rate
	Catheter malfunction	<ul style="list-style-type: none"> • Check access
High Filter Pressure Drop	Filter clotting	<ul style="list-style-type: none"> • Plan circuit replacement • Increase blood flow rate • Increase pre-filter replacement rate
	Insufficient anticoagulation	<ul style="list-style-type: none"> • Optimize anticoagulation

Ejaz AA et al. Nurs Crit Care. 2007 Mar-Apr;12(2):81-5.

Michel T et al. Curr Opin Crit Care. 2018 [Epub ahead of print]

† Adapted from Ejaz et al. (2007) and Michel et al. (2018)

Key Takeaways

Key Takeaways

- Circuit pressure trends may be used to monitor circuit function and guide interventions to prolong circuit patency and preserve treatment continuity
- Variations in CRRT modalities and prescriptions make establishing absolute pressure guidelines difficult
- Pressure trends and patterns are more clinically relevant than absolute numbers
- TMP and pressure drop provide important information regarding filter clogging and clotting
- Membrane clogging due to protein and red cell deposition is associated with a slow and continuous rise in TMP
- Circuit clotting also leads to a rise in TMP (due to loss of surface membrane of clotting fibers); this may be progressive, but can often be sudden
- A slow and continuous rise of filter pressure drop is an indicator of filter clotting (fiber occlusion)
- It is important to keep in mind that the operating points for TMP and pressure drop are recalculated during therapy

References

References

- Brain M, Winson E, Roodenburg O, McNeil J. Non anti-coagulant factors associated with filter life in continuous renal replacement therapy (CRRT): a systematic review and meta-analysis. *BMC Nephrol*. 2017 Feb 20;18(1):69.
- Dunn WJ, Sriram S. Filter lifespan in critically ill adults receiving continuous renal replacement therapy: the effect of patient and treatment-related variables. *Crit Care Resusc*. 2014 Sep;16(3):225-31.
- Ejaz AA, Komorski RM, Ellis GH, Munjal S. Extracorporeal circuit pressure profiles during continuous venovenous haemofiltration. *Nurs Crit Care*. 2007 Mar-Apr;12(2):81-5.
- Gambro. Prismaflex Operator's Manual (for use with software versions 7.xx). Order Number G5039110. 2005-2014.
- Holt AW, Bierer P, Bersten AD, Bury LK, Vedig AE. Continuous renal replacement therapy in critically ill patients: monitoring circuit function. *Anaesth Intensive Care*. 1996 Aug;24(4):423-9.
- Joannidis M, Oudemans-van Straaten HM. Clinical review: Patency of the circuit in continuous renal replacement therapy. *Crit Care*. 2007;11(4):218.
- Kakajiwala A, Jemielita T, Hughes JZ, Windt K, Denburg M, Goldstein SL, Laskin B. Membrane pressures predict clotting of pediatric continuous renal replacement therapy circuits. *Pediatr Nephrol*. 2017 Jul;32(7):1251-1261.
- Kim IB, Fealy N, Baldwin I, Bellomo R. Premature circuit clotting due to likely mechanical failure during continuous renal replacement therapy. *Blood Purif*. 2010;30(2):79-83.
- Michel T, Ksouri H, Schneider AG. Continuous renal replacement therapy: understanding circuit hemodynamics to improve therapy adequacy. *Curr Opin Crit Care*. 2018 Sep 20. doi: 10.1097/MCC.0000000000000545. [Epub ahead of print]
- Ronco C, Bellomo R, Homel P, Brendolan A, Dan M, Piccinni P, La Greca G. Effects of different doses in continuous veno-venous haemofiltration on outcomes of acute renal failure: a prospective randomised trial. *Lancet*. 2000 Jul 1;356(9223):26-30.
- Uchino S, Fealy N, Baldwin I, Morimatsu H, Bellomo R. Continuous is not continuous: the incidence and impact of circuit "down-time" on uraemic control during continuous veno-venous haemofiltration. *Intensive Care Med*. 2003 Apr;29(4):575-578.
- Zhang L, Baldwin I, Zhu G, Tanaka A, Bellomo R. Automated electronic monitoring of circuit pressures during continuous renal replacement therapy: a technical report. *Crit Care Resusc*. 2015 Mar;17(1):51-4.
- Zhang L, Tanaka A, Zhu G, Baldwin I, Eastwood GM, Bellomo R. Patterns and Mechanisms of Artificial Kidney Failure during Continuous Renal Replacement Therapy. *Blood Purif*. 2016;41(4):254-63.

Baxter

Thank you

The PRISMAFLEX Control Unit is intended for:

- Continuous Renal Replacement Therapy (CRRT) for patients weighing 20 kilograms or more with acute renal failure and/or fluid overload.
- Therapeutic Plasma Exchange (TPE) therapy for patients weighing 20 kilograms or more with diseases where removal of plasma components is indicated.

Baxter and Prismaflex are trademarks of Baxter International Inc. or its subsidiaries.



**Backup Slides:
Studies of CRRT Pressure Trends**

Study: Monitoring Circuit Function

- Study to identify causes of circuit cessation and determine predictive monitors of circuit function
 - 5 consecutive CVVHD patients, 41 consecutive filters
- Filter clotting accounted for 63% of circuit cessation
- Filter clotting was associated with a rise in TMP and decreased circuit function (decrease in molecular clearance)
- An increase in TMP of ≥ 26 mmHg predicted imminent circuit failure due to clogging





FIGURE 2: Circuit cessation in patient due to haemofilter clotting is demonstrated: rapid terminal rise in the pressure gradient across the haemofilter (P-Grad) without change in post-filter pressure and the accompanying fall in creatinine clearance (creat. clear.) and urea ratio.

Holt AW et al. Anaesth Intensive Care. 1996 Aug;24(4):423-9.

Study: Pressure Profiles During CVVH


- Retrospective review of 91 CVVH cases
- Causes of circuit disruption
 - Filter clotting 46%
 - Catheter degradation/ malfunction 36%
 - Routine 72-hour filter change 33%
 - Deaeration chamber clotting 4%
- Distinct pressure patterns in patent and disrupted circuits
 - Pressures in patent circuits remained flat or close to initial values
 - Pressures in non-patent circuits varied depending on the cause of disruption
- Circuit changes due to filter clogging were distinguished by deviation of transmembrane and effluent pressures from baseline



Ejaz AA et al. Nurs Crit Care. 2007 Mar-Apr;12(2):81-5.

Study: Circuit Pressures during CRRT

- Continuous recording of circuit pressures in single CRRT treatment
- The filter clotted at 33 hours with a high TMP (313 mmHg)
- TMP and filter pressure (PFP) progressively increased (TMP, 104–313 mmHg; PFP, 131–185 mmHg)
- Reduction in effluent pressure (EP) (absolute $\Delta = -168$ mmHg) reflects the increasing negative pressure due to constant ultrafiltrate demand across a clogging membrane



Study: TMP and ΔP Predict Clotting

- Retrospective analysis of 26 children treated with 79 CRRT circuits
- 64.6% circuits underwent unplanned filter change
- Each 1-mmHg increase in TMP or filter pressure was independently associated with 1.5% and 1.5% higher risk of clotting, respectively ($p < 0.001$)




Fig. 2. Changes in hourly transmembrane pressures (*TMP*) (*y-axis*) over time (*x-axis*, in hours) in circuits that clotted prematurely (*left panel*) as compared to circuits that lasted 3 days (*right panel*). *Black dots* indicate when a circuit clotted. *Red line* is a lowess line indicating overall trend. In circuits clotting prematurely, the transmembrane pressures increased steadily over time. In circuits lasting 3 days, the transmembrane pressures appeared to be more constant, with hourly values mostly <200 mmHg

Kakajiwala A et al. *Pediatr Nephrol.* 2017 Jul;32(7):1251-1261.