oXiris

Emergency Use Authorization for the United States

The oXiris Set has been Authorized by the FDA for the indication to treat patients with COVID-19 infection.

The oXiris Set has been authorized by FDA under EUA200164.

The oXiris Set is Authorized only for the duration of the declaration that circumstances exist justifying the authorization of the emergency use of the oXiris set under section 564(b)(1) of the Act, 21 U.S.C. § 360bbb-3(b)(1), unless the authorization is terminated or revoked sooner.

Intended Use for Patients with COVID-19

The oXiris Set is indicated for use only with the Prismaflex control unit or with the PrisMax control unit.

It is intended to treat patients 18 years of age or older with confirmed COVID-19 admitted to the ICU with confirmed or imminent respiratory failure in need of blood purification, including use in continuous renal replacement therapy, to reduce pro-inflammatory cytokine levels, who have any one of the following conditions:

- Early acute lung injury (ALI)/ early acute respiratory distress syndrome (ARDS); or
- Severe disease, defined as:
 - dyspnea,
 - respiratory frequency ≥ 30 /min,
 - blood oxygen saturation $\leq 93\%$,
 - partial pressure of arterial oxygen to fraction of inspired oxygen ratio < 300 , and/or
 - lung infiltrates $>50\%$ within 24 to 48 hours; or
- Life-threatening disease, defined as:
 - respiratory failure,
 - septic shock, and/or
 - multiple organ dysfunction or failure

The use of oXiris is contraindicated (as mentioned in the IFU) when:

- Patients present a known allergy to heparin or have type II thrombocytopenia caused by heparin (HIT Syndrome type II)
- A drug used simultaneously with oXiris is contraindicated per its Instructions for use

Relative contraindications (individual risk/benefit to be determined by treating physician) for the use of oXiris include:

- The inability to establish vascular access to safely perform CRRT / hemoperfusion (SCUF; CVVH; CVVHD; CVVHDF)
- Severe hemodynamic instability
- Known hypersensitivity to any component of the oXiris Set

This set is intended for use in the following veno-venous therapies: SCUF; CVVH; CVVHD; CVVHDF.

All treatments administered with the oXiris Set must be prescribed by a physician. The size, weight, state of uremia, cardiac status, and general physical condition of the patient must be carefully evaluated by the prescribing physician before each treatment.

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The **oXiris** set is manufactured by **GAMBRO Industries, 7 avenue Lionel Terray, BP 126, 69883 MEYZIEU CEDEX, FRANCE.**

⚠ **Caution:** Federal law (USA) restricts this device to sale by or on the order of a physician.

DEFINITION OF EXPRESSIONS USED IN THIS MANUAL

In this document :

⚠ **“Warning”** indicates a hazardous situation which, if not avoided, could result in death or serious injury.

⚠ **“Caution”** indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

“Note” to give additional information.

SCUF: Slow Continuous UltraFiltration.

CVVH: Continuous Venous-Hemofiltration.

CVVHD: Continuous Venous-Hemodialysis.

CVVHDF: Continuous Venous-Hemodiafiltration.

Predilution: addition of replacement fluid to the blood stream upstream to the filter.

Postdilution: addition of replacement fluid to the blood stream downstream to the filter.

"Control unit" refers to the **PrismaFlex** control unit or to the **PrisMax** control unit (in countries where **PrisMax** is cleared or registered).

PRODUCT DESCRIPTION

- The **oXiris** set is a disposable, extracorporeal circuit for use with the **PrismaFlex** control unit or with the **PrisMax** control unit (in countries where **PrisMax** is cleared or registered).
- The **oXiris** set consists of a hollow fiber hemofilter/dialyzer* with an internal surface permanently grafted with highly purified heparin from porcine origin, and tubing lines; refer to the control unit operator manual drawing for details.
- The **oXiris** set enables removal of various inflammatory mediators by diffusion, convection and adsorption, including cytokines and endotoxins.
- The blood return line (blue-striped) is equipped with a Luer-lock connection near the deaeration chamber, dedicated to the connection of authorized devices and accessories described in the control unit operator's manual.
- All line connectors are compatible with the ISO 594-1 and ISO 594-2 international standards concerning conical fittings.
- The fluid pathways of the **oXiris** set are guaranteed sterile and non pyrogenic.
- The **oXiris** set is sterilized by ethylene oxide (EtO). Deaeration is such that EtO residuals comply with the standards in ISO 10993.
- Expiration date: please refer to product label.

* In this document the hemofilter/dialyzer will be referred to as "filter".

INTENDED USE / INDICATIONS

The **oXiris** set is indicated for use only with the **PrismaFlex** control unit or with the **PrisMax** control unit (in countries where **PrisMax** is cleared or registered). It is intended for patients in need of blood purification, including continuous renal replacement therapy, and in conditions where excessive endotoxin and inflammatory mediators levels exist.

This set is intended for use in the following veno-venous therapies: SCUF; CVVH; CVVHD; CVVHDF.

All treatments administered with the **oXiris** set must be prescribed by a physician. The size, weight, state of uremia, cardiac status, and general physical condition of the patient must be carefully evaluated by the prescribing physician before each treatment.

CONTRAINDICATIONS

It is contra-indicated to use the **oXiris** set where patients present a known allergy to heparin or have type II thrombocytopenia caused by heparin (HIT Syndrome type II).

Any contra-indications indicated in the respective Instructions for use of all drugs used simultaneously with this set should be taken into account.

For the following conditions a careful assessment of the individual risk/benefit ratio has to be made by the treating physician (relative contraindications):

- inability to establish vascular access,

- severe hemodynamic instability,
- known hypersensitivity to any component of the **oXiris** set.

CAUTIONS AND WARNINGS

Note: refer to the control unit user interface and operator's manual for additional cautions and warnings.

⚠ Cautions

1. Particular attention must be paid to extra corporeal blood volume with respect to patient size. Consider the sum of the **oXiris** set blood volume (refer to "Specifications") plus the blood volume of any accessory or device if used. The **oXiris** sets should be restricted to patients with a body weight greater than 30kg (66lb).
2. If the patient is not immediately connected to the **oXiris** set after priming is complete, flush the set with at least 1 000 mL priming solution (without heparin added) prior to connecting the patient. This requires use of a new bag of priming solution.
3. When not using the pre blood pump infusion line, it is recommended to clamp this line close to its connection to the access blood tubing line; this will prevent the sedimentation of blood into the pre-blood pump infusion line.
4. Using the **oXiris** set with blood flow rates lower than the recommended minimum values (see "Operating Parameters" section) may impair filter performance due to hemoconcentration, or to increased risk of coagulation.
5. Since drugs can be removed by the membrane of the filter, the dosage of associated drug treatments may need to be adjusted for patients on continuous renal replacement therapy. Monitoring of blood drug levels of relevant compounds should be performed. The removal of other water-soluble compounds (e.g. vitamins, trace elements) during therapy also requires clinical consideration.
6. The membrane is pre-heparinised, however it possesses residual capacities for the adsorption of heparin. In order to saturate the membrane it is recommended to use a heparinised rinsing solution. The life span of the set could be affected if this procedure is not respected.

⚠ Warnings

1. Carefully read these instructions for use and the control unit operator's manual before using this product.
2. The use of operating procedures other than those published by the manufacturer or the use of accessory devices not recommended by the manufacturer can result in patient injury or death.
3. Store the **oXiris** set in a dry place, between 0° C (32° F) and 30° C (86° F).
4. Do not use this set if the packaging is damaged, if the sterilization caps are missing or loose, or if any of the lines in the set are kinked.
5. To prevent contamination, this **oXiris** set must be used as soon as its packaging and sterilization caps are removed.
6. Do not try to remove the filter from the cartridge plate.
7. Use aseptic techniques when handling all blood and fluid lines in the set.
8. **oXiris** sets are compatible with the usual disinfection agents used for aseptic setup; however solvents and other chemicals, if used in contact with the product, could damage the set.
9. During priming and operation, observe closely for leakage at joints within the set, and connections to other approved accessories and bags. Leakage can cause blood loss, fluid imbalance or air embolism. If a leakage is detected at a Luer connection and cannot be stopped by tightening the connections, or if leakage occurs at any other location, replace the set.
10. Tightening Luer connections with an excessive force can damage the connectors.
11. Once priming is complete the **oXiris** set's blood circuit will still contain heparinised saline solution. Depending on the level of the patient's bleeding risk the physician must decide if an additional priming using 500 mL non heparinised saline solution is necessary.
12. Do not allow air to enter the blood compartment of the filter after priming is started. If a large amount of air enters, the set must be replaced.
13. Should acute allergic reactions (first-use syndrome) occur in patients receiving treatment, **immediately stop the treatment** and administer appropriate intervention. Pay special attention to patients receiving ACE inhibitors and/or having already shown similar allergic reactions (see "Hypersensitivity Reactions" section).
14. Use a 21-gauge or smaller needle to obtain blood/fluid samples or remove trapped air from the **oXiris** set. Use of larger needles can cause holes in the sample sites, resulting in external leak or air intake.
15. External blood leakage may not be immediately identified by monitoring equipment and could result in significant blood loss. Check the filter and all connections of the disposable tubings during treatment to minimize the risk of leakage.

16. To assure adequate filter performance, it is recommended that the set be changed every 24 hours of use. However, the set must be changed after 3 days (72 hours). Continued use beyond this limit could result in rupture of the pump segments, with risk of patient injury or death.
17. Destroy this set after single use, using aseptic technique for potentially contaminated equipment and following local regulation for disposal. Do not re-sterilize. The **oXiris** set is intended for single use only. Re-using the **oXiris** set may cause serious damage to the product resulting in patient injury or death.
18. Use only drugs compatible with plastics listed in the specifications section. Some plastics can be incompatible with drugs when in contact with solutions with pH > 10.

SPECIFICATIONS

See Tables at end of document.

SET MATERIALS

AN69 HF hollow fiber :	Acrylonitrile and sodium methallyl sulfonate copolymer + Polyethylenimine (surface treatment agent) + heparin grafted
Housing and headers :	Polycarbonate
Potting compound :	Polyurethane
Tubing material :	Plasticized polyvinyl chloride (PVC)
Cartridge :	Polyethyleneterephthalate Glycol

Note: the following information is available from the manufacturer upon request:

- information about test methods used to obtain performance characteristics,
- the number and range of particles in the effluent from the dialyzer prepared as recommended for clinical use,
- the types and amounts of residue from the sterilization process.

Note: the **oXiris** set is not made with natural rubber latex.

Note: the **oXiris** set is DEHP-free.

INSTRUCTIONS FOR USE

Note: use the set by following the detailed on-line instructions provided by the control unit. Additional information is available in the control unit operator's manual.

Note: In case no renal replacement therapy is prescribed by the physician, it is recommended to select the SCUF mode to limit the number of bags used. The resulting modality will be equivalent to hemoperfusion in case no patient fluid removal is prescribed.

Load Set

Install the set onto the control unit using the photographs on the inside cover as a guide.

Prepare and Connect Solutions

The membrane is pre-heparinised, however it possesses residual capacities for the adsorption of heparin. In order to saturate the membrane it is recommended to use a heparinised rinsing solution. The life span of the set could be affected if this procedure is not respected.

Hang bag of priming solution with added 5000 IU unfractionated heparin/liter (correctly homogenized) on priming hook.

Connect access (red)/effluent (yellow) Y-line to priming solution bag.

SPECIAL PROCEDURES IN CASE OF COMPLICATION

External Blood Leaks

Note: see Warning no. 15.

If an external blood leakage is observed, immediately stop the blood pump. Initiate corrective action by securing connections or replacing the **oXiris** set.

If necessary, administer adequate replacement solution to the patient to compensate for blood loss.

Hypersensitivity Reactions

Note: see Warning no. 13.

Should acute allergic reactions (first use syndrome) occur within the first few minutes of the treatment, it is important to react immediately by discontinuing the session and administering appropriate treatment.

Adverse reactions may occur due to the complex interaction between blood and the artificial surfaces of the entire extracorporeal circuit. These reactions may also be precipitated and/or exacerbated by other external factors involved with the individual patient's specific disease process and the treatment of renal insufficiency. Certain types of adverse reactions may occur due to operational factors associated with the treatment. Therefore, proper management of the fluid removal, electrolyte balance, anticoagulation and blood flow rate as well as monitoring of the overall treatment parameters are essential to avoid side-effects which may be associated with hemodialysis/hemofiltration therapies.

Hypersensitivity reactions have been observed during dialysis. Symptoms of a hypersensitivity reaction may be gastrointestinal, mucocutaneous, respiratory, cardiovascular or systemic in nature and range from very mild to severe. Such symptoms have been described as anaphylactic-like reactions within the first few minutes. Manifestations include nausea, malaise, weakness, a sensation of burning or heat throughout the body, profuse perspiration, respiratory distress and in some instances hypotension and cardiopulmonary arrest. Should a combination of such symptoms appear, particularly at the start of the treatment session, it is important to react immediately by discontinuing the session and administering appropriate treatment. Blood in the extracorporeal circuit must not be returned to the patient.

Extra care must be taken when treating patients who have exhibited possible hypersensitivity symptoms during previous treatments, or patients who have a history of being highly sensitive and allergic to a variety of substances. A physician must be consulted to evaluate the risk and prescribe the appropriate precautions if a possible sensitivity is suspected.

The following factors are considered essential to minimize the risk of hypersensitivity reaction and other side effects:

- Strict adherence to the set-up, priming and rinsing procedures detailed in the manufacturer's instructions for use.
- Setting up and monitoring the treatment operating parameters according to the manufacturer's recommendations specified for each type of **oXiris** set and to the patient's needs and tolerance.
- Strict adherence to all WARNINGS and CAUTIONS given by the manufacturer in the instructions for use.

Patients receiving angiotensin converting enzyme (ACE) inhibitors as medication can develop, within the first few minutes of a treatment, symptoms similar to acute allergic reactions i.e bronchospasm, edema of airways or larynx, dyspnea, angioedema, urticaria, nausea, vomiting, diarrhea, respiratory arrest, abdominal cramping, hypotension, hypovolemic shock and death.

However, for these patients, administration of antihistamines often does not alleviate the symptoms. In this case, treatment must be stopped and a more aggressive first-line therapy for an anaphylactoid reaction should be initiated immediately after the onset of symptoms.

Therefore, the prescribing physician must pay special attention to patients receiving ACE inhibitors and/or having already shown similar reactions.

WARRANTY AND LIMITATION OF LIABILITY

- a) The manufacturer warrants that the **oXiris** set has been manufactured in accordance with its specifications and in compliance with good manufacturing practices, other applicable industry standards and regulatory requirements.
If provided with the lot/serial number of the defective product, the manufacturer will, by replacement or credit, remedy manufacturing defects in the **oXiris** set becoming apparent before the expiration date.
- b) The warranty under paragraph a) above is in lieu of, and to the exclusion of, any other warranty, whether written or oral, express or implied, statutory or otherwise, and there are no warranties of merchantability or other warranties, which extend beyond those described in paragraph a) above. The remedy set out above for manufacturing defects is the sole remedy available to any person due to defects in the **oXiris** set and the manufacturer shall not be liable for any consequential or incidental loss, damage, injury or expense arising directly or indirectly from the use of the **oXiris** set, whether as a result of any defect therein or otherwise.
- c) The manufacturer shall not be liable for any misuse, improper handling, non-compliance with warnings and instructions, damage arising from events after the manufacturer's release of the **oXiris** set, failure or omission to inspect the **oXiris** set before use in order to ensure that the **oXiris** set is in proper condition, or any warranty given by independent distributors or dealers.
- d) The manufacturer is GAMBRO Industries, 7 avenue Lionel Terray, BP 126, 69883 MEYZIEU CEDEX, FRANCE.

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⁽¹⁾ Nominal values - given for indication

⁽²⁾ Typical mean values obtained from laboratory testing of post-sterilization sample lots. Results may vary depending on patient and clinical conditions. Adsorption removal rate obtained in vitro are likely to differ from in vivo results. Adsorption characteristics change with the duration of observation.

⁽³⁾ Ultrafiltration is controlled by the control unit and is independent of the ultrafiltration coefficient (KUF). Note : a TMP > 40 kPa (300 mmHg) does not allow a higher ultrafiltration

⁽⁴⁾ Removal Rate expressed at t=120min with a theoretical initial IL-10, IL-6, HMGB-1 and TNF-α respective concentration of 500 pg/mL, 1500 pg/mL, 30 ng/mL and 250 pg/mL.




⁽⁵⁾ Removal Rate expressed at t=120min with an initial LPS concentration after stabilization of 50±10 EU/mL.




Prismatherm

SP420 line connexion

INSTALL SET DIRECTLY ONTO THE CONTROL UNIT

Patented ready to use PrismaFlex set packaging.